

PRODUCT REVISION ANNOUNCEMENT

RELEASE DATE: 25/06/2020 - Document DSP0620

This announcement contains revision of Dispensers Vol2, in Health Protection Equipment Index of MT Displays 2020 Range.

Dear MT Displays Partners, Distributors and Other Interested Partners,

Effective from today, old technology Dispensers Vol2 have been revised. Through this specification development process, MT Displays ensures quality of dispenser. As a key partner, MT Displays invites you to participate in following below information to ensure continuity of Dispensers Vol2.

The tentative MT Displays improvement revision schedule includes the following Milestones:	Dates
Release specification announcement letter	Week 26, 2020
Coordinate with distributors for pending orders	Ongoing
Release new technical specifications	June, 2020

Steps of Revision:

- The Dispenser Vol2 has been designed to improve its functionality.
- Dispenser holder is made of shutter-proof ABS to allow inserting into the profile or post directly.
- Spring loaded locking system prevents foreign access to the bottle.
- Lever stays open for easy change of dispenser.
- Pump is placed into the nest under the lever, to avoid pump to move around.
- Pump nozzle faces downwards allowing alcohol/gel into hand, without wasting sanitizing gel/alcohol or splashing around.
- Drip tray inserts into the holder with a basic locking system.

BOM Changes

Existing SKUs which include dispenser Vol2 have remained the same. Product trees have been revised automatically.

Distributor Participation:

The new and old dispensers are not replaceable. Product in complete or parts cannot be mixed for new and old stock. Products scheduled for shipment as of Week 26, 2020 will be shipped with the new improved dispenser Vol2.

